

# Software Installation Instructions for AlcoMate AccuCell (AL9000) on Windows 7 / Vista

1. DO NOT PLUG THE USB CABLE INTO YOUR PC YET!
2. Insert the MAX20 CD into your PC. The serial number on the CD Sleeve will be used later, after MAX20 software is installed.
3. Navigate to the CD drive on your PC and run (double-click) the file starting with "CP210x". Follow the onscreen instructions to proceed by clicking "Next" and accepting license terms.
4. Toward the end of the setup, you MUST place a check in the box to "Launch the CP210x VCP Driver Installer" and then click "Finish". Please make sure you do not skip this step.
5. Install the CP210x USB to UART Bridge Driver.
6. After successful installation, connect the USB cable to your device and then to your PC. You MUST wait for your computer to display a message that the device is either "successfully installed" or "ready to use". Please make a note of the COM port (for example, COM3 or COM4) that is assigned, if troubleshooting is necessary.
7. Go back to the CD drive on your computer and run (double-click) the "Setup" file to install MAX20.
8. After MAX20 is installed, navigate to the MAX20 program from your Start Menu (DELETE the shortcut on your desktop). RIGHT-CLICK on the MAX20 program icon in your Start Menu and select "Properties".
9. In the Properties window, select the "Compatibility" tab. Place a check in front of "Run this program in compatibility mode for:" and select "Windows XP (Service Pack 2)" from the drop-down box. Also place a check in front of "Run this program as an administrator". Click "Apply" then "OK".
10. You can now run the MAX20 program from the Start Menu. In the Registration tab, you can create ID numbers and Names for people being tested. There are 4 default entries. In the Measure tab, you can select one of these IDs, then begin testing. If you are able to click the "Start" button and your device indicates "blow", then everything is functioning properly.